


## Anorectal Balloon Expulsion Test

### Anorectal Balloon Expulsion Catheters

The anorectal balloon expulsion catheter is used to assess the patient's ability to defecate.

- ✓ Used for Balloon Expulsion Test and Biofeedback.
- ✓ Single-Use
- ✓ Specially designed for easy placement.
- ✓ Non latex 400mL volume balloon
- ✓ Measurement lines for balloon placement.
- ✓ Can be used with Solid State or water perfusion motility systems



#### Please Pay Attention!

This balloon expulsion catheter is designed for single use.

Re-sterilizing and/or attempting reuse may result in device failure and/or spread of disease. To prevent reuse, the balloon expulsion catheter should be disposed of immediately after use.

The catheter and its packaging should be checked for any mechanical damage or malfunction.

If damaged, it should not be used. Silicone-based lubricant should not be used with the balloon expulsion catheter. The shelf life and duration of use of the Anorectal Balloon Expulsion catheter is two years.

Standard Single-Use Anorectal Balloon Expulsion Catheter

REF SR1B



Made in Canada.

